

PEROXAN BD-Paste 50 SI

Crosslinking Diacyl peroxides

Description: Di-(2,4-dichlorobenzoyl)-peroxide

50%, Paste in silicone oil

PEROXAN BD-Paste 50 SI is used for the crosslinking of silicone rubber.

Molecular weight: 380,0 CAS No.: 133-14-2

Technical data: Appearance: white paste

Peroxide assay: ca. 50% Active oxygen assay: ca. 2,10%

Density at 20°C: 1,18 g/cm³

Half life time: in a silicon rubber compound

Storage: Maximum storage temperature (T_{s max}): 30°C

Minimum storage temperature ($T_{s min}$): 5°C Storage stability as from date of delivery: 6 months

Organic Peroxides are more or less stable products but will decompose under the influence of heat. To minimize a loss of quality during storage, it is important that the recommended maximum storage temperature is not exceeded. If a minimum storage temperature is given, an undesirable process such as a solidification or phase separation, is known to occur below this temperature.

Thermal stability: SADT: 60°C

Emergency temperature (T_{em}): 55°C

The SADT is the lowest temperature at which a self accelerating decomposition may occur. The emergency temperature is derived from the SADT. It is the temperature at which emergency actions have to be taken.



PEROXAN BD-Paste 50 SI

Crosslinking Diacyl peroxides

CL/102500/01.09/2004

Page 2 of 2

Application: PEROXAN BD-Paste 50 SI is recommended for the crosslinking of silicone rubber.

With PEROXAN BD-Paste 50 SI silicone rubber compounds can be cured without external

pressure (hot air and/or IR vulcanization).

PEROXAN BD-Paste 50 Si can be incorporated easily into a silicone rubber compound

on a 2-roll mill.

Dosing (silicone rubber): 1,1 to 2,3 phr

Safe processing temperature (t_2): 75°C Typical crosslinking temperature (t_{90}): 90°C

The safe processing temperature t_2 is defined as the temperature, at which the scorch time is longer than 20 minutes. The typical crosslinking temperature t_{90} is defined as the temperature at which 90% of the crosslinks in the compound are formed within about 12 minutes.

Packaging: 20 kg Pail

Major decomposition Carbon dioxide, 1,3-Dichlorobenzene, 2,4-Dichlorobenzoic acid, Traces of

products: 2,2',4,4' tetrachlorobiphenyl

Safety and handling: Please refer to the material safety data sheet (MSDS) for information concerning safe storage, use

and handling of PEROXAN BD-Paste 50 SI. This information should be thoroughly reviewed prior to acceptance of this product. The MSDS is available for downloading at www.pergan.com or through

Telephone: +49 2871 9902-0

Telefax: +49 2871 990250

www.pergan.com

contacting Pergan directly.

The information presented herein is true and accurate and to the best of our knowledge, but without any guarantee. Since the conditions of use are beyond our control we disclaim any liability, including for patent infringement, incurred in connection with the use of these products, data or suggestions.